


MIX & MATCH

**Buy 10
Spoos
Get 1
Spool
FREE***

RMO's® Exclusive Energy Chain™

RMO® is excited to introduce the new Adjoined spacing option for Energy Chain™. The Adjoined spacing option is positioned between Closed and Reduced sizes. Energy Chain™ is the undisputed leader among all brands of elastic chain. University and clinical researchers have shown that RMO® Energy Chain™ consistently outperforms all other elastic chains.



Energy Chain™ retains a higher percentage of original force than Industry Leaders†


■ Initial Force (In grams)
 ■ Remaining Force at 5 weeks (In grams)

† Studies show competitor chain materials lose 50% of their initial force after 72 hours or less, while Energy Chain™ retains more than 80% of its original force after 5 weeks.



NEW!

Ring Size	Small	Small	Small	Small	Large
Space Between Rings	Closed	Adjoined	Reduced	Medium	Large
Space Between Centers	2.8 mm	3.0 mm	3.8 mm	4.0 mm	5.1 mm
Clear	J00120	J00202	J00269	J00121	J00122
Gray	J00166	J00203	J00270	J00167	J00168
Blue	J00223	J00205	J00271	J00224	J00281
Green	J00227	J00206	J00272	J00228	J00282
Red	J00231	J00204	J00273	J00232	J00280
Purple	J00235	J00207	J00274	J00236	J00283
Teal	J00136	J00214	J00137	J00138	J00290
Pearl Blue	J00139	J00215	J00140	J00141	J00291
Orange	J00239	J00208	J00275	J00240	J00284
Porsche Red	J00261	J00209	J00276	J00262	J00285
Black	J00263	J00210	J00277	J00264	J00286
Light Pink	J00265	J00211	J00278	J00266	J00287
Yellow	J00267	J00212	J00279	J00268	J00288
Ivory	J00142	J00213	J00143	J00144	J00289
Pink	J00145	J00216	J00146	J00147	J00292
Clear	J00293	J00295	J00297	J00299	J00301
Gray	J00294	J00296	J00298	J00300	J00302



* All Buy 10 Get 1 Free and Mix and Match products must be of equivalent products to qualify. ** Prices are subject to change at any time without notice.

Visit RMOOrtho.com for complete catalog of products.