

FULLY COMPATIBLE ACTIVE TORQUE CONTROL

# Anterior Root Torquing Spring



## Universal compatibility — no adjustments necessary

Pre-formed torquing auxiliary  
engages the arch wire

Provides lingual or labial  
torque to the maxillary or  
mandibular anterior teeth

Universal design compatible  
with any bracket system

Requires no forming or  
bending

Torque control remains active  
throughout treatment with no  
necessary adjustments


Part numbers:

A00431 – Maxillary – Small – 10 Pkg.

A00432 – Maxillary – Medium – 10 Pkg.

A00433 – Maxillary – Large – 10 Pkg.

A00440 – Mandibular – Universal – 10 Pkg.



**RMO**<sup>®</sup>

**RMOrtho.com**  
customerservice@rmortho.com  
2165 Earlywood Dr. | Franklin, IN 46131

Toll Free: 800.525.6375  
Fax: 303.592.8200  
Local: 303.592.8209



Scan here to order your  
favorite RMO<sup>®</sup> products!  
Order online at **RMOrtho.com**