

# RMO® Direct Bond Attachments

Check out our **NEW** Bondable Products!



### Eruption Appliance

- Plated stainless steel chain and base
- 80 gauge mesh pad
- Adjustable 2cm chain

Silver, Pack of 5 .....A00310

Gold, Pack of 5 ..... A00311



### Direct Bond, Bite Guide

- One-piece design
- Efficiently open a deep bite
- Achieve mandible advancement

Pack of 10 .....A00317



### Tongue Trainer

- Curved base
- Help control tongue thrusting habit

Pack of 10 .....A00318



### Direct Bond Button, Round Ceramic

- Great for use with aligners or aesthetic bracket systems

Pack of 10 .....A00321



### Molded Plastic Lingual Button

Pack of 10 ..... A03902



### Direct Bond Eyelet, Rectangular

- One-piece curved base

Pack of 10 .....A00314



### Direct Bond Coil Eyelet, Round

- Round base

Pack of 10 .....A00315



### Traction Hook, Round

- 80-gauge mesh pad

Pack of 10 .....A00316



### Direct Bond Lingual Cleat

- Stainless steel
- Low Profile
- Easily Adjustable

Pack of 10 ..... A00407



### Bondable Lingual Cleat

- Low profile design reduces irritation

Pack of 10 ..... A07911



### Direct Bond Lug/Cleat (Debnam)

- Mesial-distal elastic hook
- Hook is bent out from base surface to engage elastics

Pack of 10 .....A07910



### Direct Bond Eyelet

- Accepts wire and small diameter elastic ligatures

Pack of 10 ..... A07914

### Bondable Buttons

Stainless steel low-profile buttons on 80 mil mesh base.

- Round or rectangular base
- Flat or curved contour bonding base
- Smooth finish with rounded edges for patient comfort



### Round Curved Base

Pack of 10 ...A07916



### Rectangle Curved Base

Pack of 10 ...A07918



### Round Flat Base

Pack of 10 ... A07917



### Rectangle Flat Base

Pack of 10 ...A07919

M0824006 REV. A



RMOOrtho.com  
CustomerService@RMOOrtho.com  
2165 Earlywood Dr. | Franklin, IN 46131

Toll Free: 800.525.6375  
Fax: 303.592.8200  
Local: 303.592.8209

Scan here to order your favorite RMO® products!  
Order at [RMOOrtho.com](http://RMOOrtho.com)

