

Apex Series™ Aligner Pliers




Actual Size

Surgical-Grade
Stainless Steel


Ergonomic Handles

Quickly add features
to any aligner for extra
control over corrections




Hole Punch
For relief around bonded button or
bracket.

T00571




Tear Drop Hole
To add Elastics Hook.

T00572




Horizontal Indent
Indent to add root torque.

T00573




Vertical Indent
Create mesial or distal indent for
rotation over correction.

T00574



Eraser Plier
Returns aligner surface to original
shape post adjustment.

T00576



Retention Dimple
Creates round dimple to add
retention as desired.

T00575

Apex™ Series is trademarks of OrthoAmerica Holdings, LLC.
© 2025 OrthoAmerica Holdings, LLC.

M0824007 REV. A

RMO®

RMOrtho.com
CustomerService@RMOrtho.com
2165 Earlywood Dr. | Franklin, IN 46131

Toll Free: 800.525.6375
Fax: 303.592.8200
Local: 303.592.8209

Scan here to order your
favorite RMO® products!
Order at **RMOrtho.com**

