



# FLI<sup>®</sup> Clear M Brackets


## Aesthetics and Performance Personified

### The Clear Choice in Patient Aesthetics

- FLI Clear M monocrystalline brackets offer exceptional purity and integral strength
- Aesthetically attractive and patient friendly with contoured edges and low profile height
- These sapphire grade brackets are outstanding alternatives to conventional metal brackets without compromising performance
- All FLI Clear M brackets incorporate alumina oxide powder coated bonding bases creating reliably strong mechanical adhesion, ensuring endurance yet allowing for safe and predictable debonding



PRESCRIPTIONS: Roth, MBT<sup>™</sup>  
SLOT SIZE: .018 & .022



#### Colored Orientation Guides

Each FLI Clear M Bracket features removable colored orientation guides allowing for accurate placement on the long axis and occlusal plane.


## Roth Prescription – Maxillary Brackets (Pack of 5)

Tooth	Torque	Angulation	Offset	Hook	.018		.022	
					Right	Left	Right	Left
Central	11°T	5°	0°		Q27200	Q27201	Q27300	Q27301
Lateral	8°T	9°	0°		Q27202	Q27203	Q27302	Q27303
Cuspid	-2°T	8°	0°	Distal	Q27204	Q27205	Q27304	Q27305
1st/2nd Bicuspid	-7°T	0°	0°		Q27206	Q27207	Q27306	Q27307
	-7°T	0°	0°	Distal	Q27208	Q27209	Q27308	Q27309



## Roth Prescription – Mandibular Brackets (Pack of 5)

Tooth	Torque	Angulation	Offset	Hook	.018		.022	
					Right	Left	Right	Left
Central/Lateral	0°T	0°	0°		Q27210		Q27310	
Cuspid	-11°T	7°	0°	Distal	Q27211	Q27212	Q27311	Q27312
1st Bicuspid	-17°T	0°	0°		Q27213	Q27214	Q27313	Q27314
	-17°T	0°	0°	Distal	Q27215	Q27216	Q27315	Q27316
2nd Bicuspid	-22°T	0°	0°		Q27217	Q27218	Q27317	Q27318
	-22°T	0°	0°	Distal	Q27219	Q27220	Q27319	Q27320



## Roth Prescription – Kits 5x5

Description	.018	.022
FLI CLEAR M ROTH Upper/Lower 5-5 Hooks on 3	KQ27201	KQ27301
FLI CLEAR M ROTH Upper/Lower 5-5 Hooks on 3,4,5	KQ27202	KQ27302




## MBT Prescription – Maxillary Brackets (Pack of 5)

Tooth	Torque	Angulation	Offset	Hook	.018		.022	
					Right	Left	Right	Left
Central	17°T	4°	0°		Q27400	Q27401	Q27500	Q27501
Lateral	10°T	8°	0°		Q27402	Q27403	Q27502	Q27503
Cuspid	0°T	8°	0°	Distal	Q27404	Q27405	Q27504	Q27505
1st/2nd Bicuspid	-7°T	0°	0°		Q27406	Q27407	Q27506	Q27507
	-7°T	0°	0°	Distal	Q27408	Q27409	Q27508	Q27509



## MBT Prescription – Mandibular Brackets (Pack of 5)

Tooth	Torque	Angulation	Offset	Hook	.018		.022	
					Right	Left	Right	Left
Central/Lateral	-6°T	0°	0°		Q27410		Q27510	
Cuspid	0°T	3°	0°	Distal	Q27411	Q27412	Q27511	Q27512
1st Bicuspid	-12°T	0°	0°		Q27413	Q27414	Q27513	Q27514
	-12°T	0°	0°	Distal	Q27415	Q27416	Q27515	Q27516
2nd Bicuspid	-17°T	0°	0°		Q27417	Q27418	Q27517	Q27518
	-17°T	0°	0°	Distal	Q27419	Q27420	Q27519	Q27520



## MBT Prescription – Kits 5x5

Description	.018	.022
FLI CLEAR M MBT Upper/Lower 5-5 Hooks on 3	KQ27401	KQ27501
FLI CLEAR M MBT Upper/Lower 5-5 Hooks on 3,4,5	KQ27402	KQ27502



\* All Buy 10 Get 1 Free and Mix and Match products must be of equivalent products to qualify.

\*\* Prices are subject to change at any time without notice.

© 2025 OrthoAmerica Holdings, LLC.

M1024017 REV. A



**RMOOrtho.com**  
CustomerService@RMOOrtho.com  
2165 Earlywood Dr. | Franklin, IN 46131

**Toll Free: 800.525.6375**  
Fax: 303.592.8200  
Local: 303.592.8209

Scan here to order your favorite RMO® products!  
Order at [RMOOrtho.com](http://RMOOrtho.com)

