

# FlashMax P4 Ortho Pro

Ortho Pro LED Curing Light Pen

## The Most Powerful Curing Light on the Market\*



True 3-Second Cure for Maximum Efficiency



Highest output intensity of LED curing lights



Superior Battery Life

- Fully recharged in 2.4 hours
- 200+ cures on a single charge
- Rechargeable lithium-ion phosphate cell

### FlashMax P4 Design Updates

- Improved Strength
- Better Grip
- Redesigned Charging Station



\* Rated by Dental Advisor Magazine as "The most powerful curing light on the market."

## FEATURES AND BENEFITS

### Adjustable tip

Light tip can adjust  $\pm 15^\circ$  for better access to distal regions of the mouth; 8mm light tip is ideal for use with fixed appliances

### Pen grip

Economical unit provides less heat in handling without sacrificing power and efficiency

### Bilateral identical control buttons

Treat the upper and lower arches with maximum efficiency using the bilateral control buttons. Identical control buttons on either side allow you to rotate the light to achieve the best angle for maximum bond strength.

### Output

- Peak Output Intensity of 6,000 mW/cm<sup>2</sup>
- 1s/3s Curing Modes

### Cordless

Provides maximum flexibility in the office; go from chair-to-chair or office-to-office

The FlashMax P4 has three programs:

	1s / 2s	3s / 5s
<b>Manual Use</b>		
<b>Double Tap</b>	2x 1s / 2s with 0.5s pause	2x 3s / 5s with 0.5s pause
<b>Rapid Automated Function</b>	1s / 2s with 0.5s pause cont.	3s / 5s with 0.5s pause cont.

### Power supply

Input Power 100-240VAC, 50/60Hz  
Output power 5VDC, 1.0A

### Battery

Nominal voltage 3.2 V  
Nominal capacity 1200 mA  
Battery type LiFePO<sub>4</sub>  
Recharging time 2.5 hours



Continuous curing mode, 3 seconds or 5 seconds with 0.5 second pause, is ideal for direct or indirect bonding.

Weight 120g  
Dimensions L 225 mm, D 25 mm

## FlashMax P4 Light Pen & Accessories

	Part No.	Item	Quantity
	J02130	FlashMax P4 Ortho Pro	1
<b>FlashMax P4 Light Pen Covers</b> are disposable and available in bags of 10 pieces.	J02121	FlashMax P4 Light Pen 4mm Tips	50
	J02122	FlashMax P4 Light Pen 8mm Tips	50
	J02123	FlashMax P4 Light Pen Covers	10
	J02124	FlashMax P4 Light Pen Shield	1
4mm	J02125	FlashMax P4 Light Pen Disposable Light Shields	50
8mm	J02126	Disposable Light Sleeves	250

M1024018 REV. A

© 2025 OrthoAmerica Holdings, LLC.



RMOOrtho.com  
CustomerService@RMOOrtho.com  
2165 Earlywood Dr. | Franklin, IN 46131

Toll Free: 800.525.6375  
Fax: 303.592.8200  
Local: 303.592.8209

Scan here to order your favorite RMO® products!  
Order at [RMOOrtho.com](https://www.RMOOrtho.com)

