

Key-Ties™ Elastomeric Ligatures




Quality & Performance Made in the USA

- RMO self-manufactures durable, high-quality elastomerics and utilizes precise dimensional engineering to ensure all force expectations are consistently met.
- Ideal ligature engagement forces for easy placement onto bracket tie wings.
- Excellent elastic recovery and long-term retention to achieve ultimate wire engagement into bracket wire slot.



Key-Ties™

- Designed with ideal force when ligating today's most popular twin bracket designs.
- .120" Outside Diameter (O.D.)
- .050" Inside Diameter (I.D.)
- 10 ties per stick, 100 sticks per pack (PK1,000)



Color	Key-Ties™
Clear	J10336
White	J10351
Pearl	J10344
Pearl Blue	J10345
Silver	J10348
Gray	J10337
Black	J10331
Burgundy	J10335
Porsche Red	J10346
Orange	J10343
Yellow	J10352
Light Green	J10340
Bright Green	J10333
Green	J10338
Aqua Green	J10330
Teal	J10350
Sky Blue	J10349
Blue	J10332
Navy	J10342
Purple	J10347
Lavender	J10339
Light Pink	J10341
Bright Pink	J10334

