

Announcing RMO's Exclusive

Zero Corrector™

in collaboration with Dr. Peri Colino




Introducing a Clinically Efficient and Versatile Class II and Class III Corrector Compatible with Both Aligner and Fixed Appliances

Strategically Offset Bar Locates Forces Closer to Center of Resistance for Improved Bodily Movement and Reduced Tipping of Canines

Contoured Comfort Hook for Easy Elastic Engagement

Integrated Hook
Expands Treatment Options such as anterior space closure



Thoughtful Design Enhancements Ensure Optimal Performance, Ease of Placement and Patient Comfort

Rigid Yet Malleable Low Profile Arm allows customizations

Articulating Molar Pad enhances molar rotating and uprighting while preventing over rotation

Anatomical Cuspid Pad ensures ideal fit for reliable bond strength

Butterfly Shape Molar Pad allows for precise positioning, easy removal of composite flash to avoid gingival irritation



Zero Corrector™

In Collaboration with Dr. Peri Colino



Case courtesy of Dr. Peri Colino



6 Accurate Sizes Patient Sets – 2 Piece

Size*	UL/LR, UR/LL
18mm	Z02018
20mm	Z02020
23mm	Z02023
25mm	Z02025
27mm	Z02027
29mm	Z02029



Individual 1 Piece Per Package

Size*	UL/LR	UR/LL
18mm	Z01118	Z01218
20mm	Z01120	Z01220
23mm	Z01123	Z01223
25mm	Z01125	Z01225
27mm	Z01127	Z01227
29mm	Z01129	Z01229



Intro Kit – 6 Patient Sets

1 set of each size KZ1000
 Contains 1 set each: 18mm, 20mm, 23mm, 25mm, 27mm, 29mm
 2x Selection Guides

***Note:** Sizing is exact from long axis cuspid to first molar.

5x Scale Model of Zero™ Corrector for Patient Consultation

7 1/16" (179.4 mm) x 1 7/16" (36.5 mm)

REF Number

P00050



Typodont Zero™ Corrector

REF Number

P00031



Zero is a trademark of OrthoAmerica Holdings, LLC. © 2025 OrthoAmerica Holdings, LLC.



RMOOrtho.com
 CustomerService@RMOOrtho.com
 2165 Earlywood Dr. | Franklin, IN 46131

Toll Free: 800.525.6375
 Local: 303.592.8200

Scan here to order your favorite RMO® products!
 Order at RMOOrtho.com

