

# RMO® Direct Bond Attachments



## Eruption Appliance

- Plated stainless steel chain ..... and base
- 80 gauge mesh pad
- Adjustable 2cm chain

Silver, Pack of 5 ..... A00310  
Gold, Pack of 5 ..... A00311



## Direct Bond, Bite Guide

- One-piece design
- Efficiently open a deep bite
- Achieve mandible advancement

Pack of 10 ..... A00317



## Tongue Trainer

- Curved base
- Help control tongue thrusting habit

Pack of 10 ..... A00318



## Bondable Power Arm

- Versatile power arm bonds directly to tooth for attachment of elastics
- 80 Gauge mesh pad
- Great for use with aligners

Left, Pack of 10 ..... A07926

Right, Pack of 10 ..... A07927



## Direct Bond Lingual Cleat

- Stainless steel
- Low Profile
- Easily Adjustable

Pack of 10 ..... A00407



## Bondable Lingual Cleat

- Low profile design reduces irritation

Pack of 10 ..... A07911



## Direct Bond Button, Round Ceramic

- Great for use with aligners or aesthetic bracket systems

Pack of 10 ..... A00321



## Molded Plastic Lingual Button

Pack of 10 ..... A03902



## Ceramic Hook

Pack of 10 ..... A00322



## Perf Button

Pack of 10 ..... A00319



## Metal Hook

Pack of 10 ..... A00320



## Direct Bond Eyelet, Rectangular

- One-piece curved base

Pack of 10 ..... A00314



## Direct Bond Coil Eyelet, Round

- Round base

Pack of 10 ..... A00315



## Traction Hook, Round

- 80-gauge mesh pad

Pack of 10 ..... A00316



## Direct Bond Lug/Cleat (Debnam)

- Mesial-distal elastic hook
- Hook is bent out from base surface to engage elastics

Pack of 10 ..... A07910



## Direct Bond Eyelet

- Accepts wire and small diameter elastic ligatures

Pack of 10 ..... A07914

## Bondable Buttons

Stainless steel low-profile buttons on 80 mil mesh base.

- Round or rectangular base
- Flat or curved contour bonding base
- Smooth finish with rounded edges for patient comfort



Round Curved Base  
Pack of 10 ..... A07916



Rectangle Curved Base  
Pack of 10 ..... A07918

