


Elastomeric Ligatures




Quality & Performance Made in the USA

- RMO self-manufactures durable, high-quality elastomerics and utilizes precise dimensional engineering to ensure all force expectations are consistently met.
- Ideal ligature engagement forces for easy placement onto bracket tie wings. Two popular force designs:
 - .120" Ideal Force
 - .115" Optimal Firm Force
- Excellent elastic recovery and long-term retention to achieve ultimate wire engagement into bracket wire slot.




Key-Ties™

- Designed with ideal force when ligating today's most popular twin bracket designs.
- .120" Outside Diameter (O.D.)
- .050" Inside Diameter (I.D.)
- 10 ties per stick, 100 sticks per pack (PK1,000)



Opti-Ties™

- Designed for use with smaller brackets with an optimal firm feel.
- .115" Outside Diameter (O.D.)
- .045" Inside Diameter (I.D.)
- 24 ties per stick, 40 sticks per pack (PK960)



Color	Key-Ties™	Opti-Ties™
Clear	J10336	J24336
White	J10351	J24351
Pearl Ivory	J10344	J24344
Pearl Blue	J10345	J24345
Silver	J10348	J24348
Gray	J10337	J24337
Black	J10331	J24331
Burgundy	J10335	J24335
Red	J10346	J24346
Orange	J10343	J24343
Yellow	J10352	J24352
Light Green	J10340	J24340
Bright Green	J10333	J24333
Green	J10338	J24338
Aqua	J10330	J24330
Teal	J10350	J24350
Sky Blue	J10349	J24349
Blue	J10332	J24332
Navy	J10342	J24342
Purple	J10347	J24347
Lavender	J10339	J24339
Pink	J10341	J24341
Bright Pink	J10334	J24334

