

# FlashMax P4 Ortho Pro

Ortho Pro LED Curing Light Pen

## The Most Powerful Curing Light on the Market\*



True 3-Second Cure for Maximum Efficiency



Highest output intensity of LED curing lights



Superior Battery Life

- Fully recharged in 2.4 hours
- 200+ cures on a single charge
- Rechargeable lithium-ion phosphate cell


### FlashMax P4 Design Updates

- Improved Strength
- Better Grip
- Redesigned Charging Station



\* Rated by Dental Advisor Magazine as "The most powerful curing light on the market."

## FEATURES AND BENEFITS



### Output

- Peak Output Intensity of 6,000 mW/cm<sup>2</sup>
- 1s/3s Curing Modes

The FlashMax P4 has three programs:

	1s / 2s	3s / 5s
<b>Manual Use</b>	[1 bar]	[1 bar]
<b>Double Tap</b>	[2 bars]	[2 bars]
<b>Rapid Automated Function</b>	[4 bars]	[4 bars]
	1s / 2s with 0.5s pause cont.	3s / 5s with 0.5s pause cont.

Continuous curing mode, 3 seconds or 5 seconds with 0.5 second pause, is ideal for direct or indirect bonding.

### Cordless

Provides maximum flexibility in the office; go from chair-to-chair or office-to-office

### Power supply

Input Power 100-240VAC, 50/60Hz  
Output power 5VDC, 1.0A



### Battery

Nominal voltage 3.2 V  
Nominal capacity 1200 mA  
Battery type LiFePO4  
Recharging time 2.5 hours



Weight 120g  
Dimensions L 225 mm, D 25 mm

## FlashMax P4 Light Pen & Accessories

	Part No.	Item	Quantity
	J02130	FlashMax P4 Ortho Pro	1
<b>FlashMax P4 Light Pen Covers</b> are disposable and available in bags of 10 pieces.	J02121	FlashMax P4 Light Pen 4mm Tips	50
	J02122	FlashMax P4 Light Pen 8mm Tips	50
	J02123	FlashMax P4 Light Pen Covers	10
	J02124	FlashMax P4 Light Pen Shield	1
4mm	J02125	FlashMax P4 Light Pen Disposable Light Shields	50
8mm	J02126	Disposable Light Sleeves	250

